

Simultaneous Measurement of Cell Color and Vacuolar pH

< Measurement of Absorption Spectrum >

1. Pour the protoplast suspension onto a x20 poly-L-lysine pre-coated glass-bottom petri dish (35 mm ϕ) and allow to stand for 1 minute at room temperature.
2. After removing the extra suspension buffer, very carefully pour 2 ml of suspension buffer into the dish.
3. Set the petri dish on a inverted microscope and measure the absorption spectrum with a micro-spectrophotometer (MCPD-7000; Photal).

< Measurement of Vacuolar pH: Preparation of Microelectrodes >

1. Make double-barreled pipettes (1B100F-4; WPI, washed with acetone and water) with vertical pullers. Cut one of the pipettes 1.5 cm from the tip.
2. Heat the pipettes at 200°C. After 1 hour, dip the blunt end of the longer pipette in 0.1% (v/v) trimethylchlorosilane/chloroform for 1 second, then heat again for one more hour.
3. Cool down the pipettes and apply melted wax to bind them together.
4. Add 0.5 μ l of proton-ionophore solution (5% v/v Hydrogen Ionophore II-Cocktail A, 0.05% w/v nitrocellulose in THF) to the longer pipette.
5. Store the pipettes in a dry box over night.
6. Add 0.2 μ l of Hydrogen Ionophore II-Cocktail A to the longer pipette, and allow to stand at vacuum condition for 5 minutes.
7. Add 0.5 M KCl to the shorter pipette.
8. Add 0.1 M MES-Tris, 0.5 M KCl, pH 6.0 to the longer pipette.

Tips • The diameter of the tip must be optimized according to the plant cells.

- All the liquid must be back-filled to the pipettes.
- The pipettes must be free of air bubbles.

< Measurement of Vacuolar pH: Calibration and pH Measurement >

1. Insert the longer pipette into the holder, and connect the shorter pipette via salt bridge.
2. Calibrate the pH electrode using standard buffers. Measure the potential difference with a electrometer (FD223; WPI), and collect data with Chart for windows v. 5.0.1 with an A/D converter (PowerLab/4sp; ADInstruments).
3. Insert the microelectrode into the cell-color recorded protoplast using a micromanipulator and record its electric potential.
4. Recalibrate the microelectrode.
5. Calculate pH_v using the latter calibration data.

Tips • Use only electrodes that have responded well'.

- Adopt data which drift was less than constant numerical value.

< Reference >

Yoshida K, Toyama-Kato Y, Kameda K, and Kondo T (2003) *Plant Cell Physiol* 44:262-268.
Yoshida K, Osanai M, Kondo T (2003) *Phytochemistry* 63:721-726.

