

## Preparation of Vacuole from Arabidopsis Leaves

### Solutions

**Medium V**; 25 mM Tris-MES (pH 7.0) and 0.45 M mannitol

**2.5% Percoll in Medium P**; 2.5% (v/v) Percoll, 25 mM Tris-MES (pH 5.5) and 0.45 M mannitol

**20% Percoll in Medium P**; 20% (v/v) Percoll, 25 mM Tris-MES (pH 5.5) and 0.45 M mannitol

**10% Percoll in Medium V**; 10% (v/v) Percoll, 25 mM Tris-MES (pH 7.0) and 0.45 M mannitol

**0.2 M K<sub>2</sub>HPO<sub>4</sub>** (Prepare before use)


**Enzyme solution**; 1% (w/v) cellulase Onozuka-RS (Yakult, Tokyo, Japan), 0.1 % (w/v) pectolyase Y-23 (Seishin, Tokyo, Japan) in **2.5% Percoll in Medium P**. (Centrifuge before use to remove debris)

### Procedure

Step1. Isolation of protoplasts

- Collect leaves in 50 ml beaker.
- Add 30 ml **Enzyme solution**.
- Vacuum infiltration.
- Gently shake for 40 min at room temperature.
- Filtrate the protoplasts with Miracloth.
- Collect the protoplasts in 15 ml test tubes.
- Add 0.5 ~ 1 ml of **20% Percoll in Medium P** to the bottom of the tubes.
- Centrifuge at 1,000 x *g* for 5 min.
- Collect heavy-protoplast fraction (Fig. 1).


Fig. 1



## Step 2. Isolation of vacuoles from protoplasts

- Transfer 1 ml of heavy-protoplasts to new 15 ml test tube.
- Add 4 ml of **0.2 M K<sub>2</sub>HPO<sub>4</sub>**.
- Seal with Parafilm and gently mix by inverting the tube five times.
- Wait **1 min 45 sec**.
- Add 4 ml of **10% Percoll in Medium V**.
- Seal with Parafilm and gently mix by inverting the tube five times.
- Remove Parafilm and add 1~1.5 ml of **Medium V** to the top of mixture.
- Centrifuge at 800 x g for 5 min (slow acceleration, slow deceleration).
- Collect vacuole fraction (ca. 40  $\mu$ l, Fig. 2).

Fig. 2



## Notes

We can recover 5 ~ 10% vacuoles of total protoplasts with 85 ~ 95% purity.